



SFC3AF

The SFC3AF is easily recognizable by its classic design. Its large gripping area makes it an easy to handle microscope throughout the classroom and even outdoors. This series is designed with the elementary school and hobby users in mind.

Model	SFC-3AF
Optical System	Finite Optical System
Observation tube	Monocular
Inclination	Adjustable inclination
Eyepiece(s)	WF10X/14mm with pointer, lockable
Nosepiece	Triple revolving nosepiece with click stops
Objective classification	Achromatic, DIN
Objective information	4X/0.10, 10X/0.25, 40X/0.65/S
Stage	Plain stage with sample clips
Stage size	88X82mm
Condenser	Rotary disc diaphragm
Focus mechanism	Coarse and fine focusing system with tension adjustment
Illumination	Plano-concave mirror and fork
Accessories included	Dust cover
Length X Width X Height	158X111X290
Weight KG	1,3KG